

Notes on Correction factors in FluorEssence

From phone conversation with Jim Mattheis (Jobin Yvon). 800 438 7739 x 122,
jim.mattheis@jobinyvon.com

Correction Factors:

R1 compensates for change or drift in lamp intensity.

- When taking an Excitation scan, you must correct with R1. ($S1/R1$)
- When taking an Emission scan, it's good to correct if experiment were being conducted over a long period of time (hours). Instrument is pretty stable—expect 1.5% change over ~5h period.

Correction files:

- Xcorr: Don't need to use this--ever. Wavelength range is narrower than R1 and the UV range is not that good.
- MCorr: Valid over ~290-850nm range. Corrects for instrument response not dark counts. Useful if comparing spectra. Not something you would use if you are measuring at a fixed wavelength.

Dark offset doesn't change significantly over time.

Good idea to record. Do this by going to Real Time Control and measuring signal with the shutter closed. Usually ~1000 cps. If you subtract dark offset from blank, also must subtract dark offset from sample scan.

Jim suggested always collecting raw data and correcting if necessary after the fact.