

# DNA Extraction Bench Protocol

This is an abbreviated protocol. Make sure to label all tubes and change pipette tips between samples.

## Sample Preparation

- Remove the abdomen of the arthropod and cut off a small piece (roughly ~2 mm, or small enough to fit in the bottom of a microcentrifuge tube). If the specimen is smaller than a grain of rice, use the entire body.
- Place the specimen in a labeled 1.5 ml microcentrifuge tube.

## Cell Lysis & DNA Precipitation


- Add 180 ul Buffer ATL and use a sterile pestle to grind the sample for 1 minute.
- Add 20 ul Proteinase K.
- Add 200 ul Buffer AL and immediately mix by vortexing for 10 seconds or pipetting up and down.
- Incubate for *at least* 15 minutes @ 56 °C. Longer incubation times (i.e., 2-3 hours) are most effective for cellular lysis.
- If arthropod debris is present, do a quick spin (~30 seconds) to pellet debris. Use a pipette to carefully transfer the supernatant to a new labeled 1.5 ml microcentrifuge tube. Discard the tube of cellular debris.
- Add 200 ul ethanol (96-100%) and mix by vortexing for 10 seconds or pipetting up and down.



Optional stopping point: store DNA in refrigerator (4 °C)

## DNA Purification

- Label a DNeasy spin column fitted with a 2.0 ml collection tube.
- Pipet the liquid containing ethanol-precipitated DNA into the DNeasy spin column.
- Centrifuge for 1 minute at  $\geq 6,000 \times g$  (8,000 rpm). Discard the flow through from the 2.0 ml collection tube.
- Add 500 ul of Buffer AW1 and centrifuge for 1 minute at  $\geq 6,000 \times g$  (8,000 rpm). Discard the flow through.
- Add 500 ul Buffer AW2 and centrifuge for 3 minutes at  $20,000 \times g$  (14,000 rpm).



## DNA Elution

- Transfer the spin columns to a labeled 1.5 ml microcentrifuge tube. Discard the 2.0 ml collection tube.
- Pipet 100 ul of Buffer AE directly onto the spin column membrane.
- Incubate at room temperature for 1 minute.
- Centrifuge at  $\geq 6,000 \times g$  or 8,000 rpm for 1 minute.
- Discard the spin column and KEEP the labeled 1.5 ml tube.
- Optional:* incubate the DNA for 1 hour @ 65 °C or overnight at room temperature.
- Store the eluted DNA frozen at -20 °C until PCR.

