
Sublimation/Rehydration of 2,5-DHA for imaging proteins

Health and Safety:

Wear nitrile gloves, laboratory coat and safety glasses. The process should be done with fume hood ventilation.

Equipment and Materials:

Equipment

Petri dish (plastic, Fisher, 08-757-100D, 100mm x 15mm); stainless steel plate (AB, 4347680); heat conductive copper tape (Electron Microscopy Sciences, 77801); Tape (Fisher, 35901R) or petri seal tape (Fisher, 50-820-970); Oven set at 37°C.

Chemicals

Trifluoroacetic acid; methanol; MilliQ water; Carnoy's fluid (6 ethanol, 3 chloroform, 1 acetic acid)

Procedure

1. Tissue preparation prior to sublimation.
Washing conditions for detection of proteins: 70% ethanol (30s), 100% ethanol (30s), Carnoy's fluid (2min), 100% ethanol (30s), H₂O (30s), 100% ethanol (30s).
2. Sublimation of 2,5-DHA:
DHA is very volatile, it can be sublimated at 105°C for 5 min to obtain ~0.15 mm/cm² at 50mm Torr.
3. Rehydration of sublimated DHA for imaging proteins (Figure 1).
 - 1) The slide with sublimated matrix is attached to a stainless-steel plate (used as a heat sink). (Note: this plate can be any type as long as it fits the size of the petri dish).
 - 2) The plate is attached to the inside of the top part of the Petri dish using a heat conductive copper tape.
 - 3) A piece of filter paper is placed in the bottom part of the Petri dish, and 500 µL of MilliQ water and 500 µL of TFA pipetted onto the paper.

- 4) The two parts of the Petri dish are reassembled so as to form a hydration chamber, sealed using a tape (Fisher, 35901R) or using petri seal tape (Fisher, 50-820-970) for the better sealing and left in the 37°C oven for 3 min.
- 5) The slide is taken out from the petri dish and is air dried.




Figure 1 Construction of the rehydration/recrystallization chamber with a Petri dish, magnet, stainless steel plate and tape.

Expected outcome:

The sign of a successful rehydration is determined by a quick color change at the moment the Petri dish is dissembled. This phenomenon is caused by the evaporation of solvent condensed on the surface of matrices.

Reference: (this is a modified procedure from following references)

1. Matrix Sublimation/Recrystallization for Imaging Proteins by Mass Spectrometry at High Spatial Resolution, Junhai Yang and Richard M. Caprioli; *Analytical Chemistry* **2011** 83 (14), 5728-5734
2. High spatial resolution imaging mass spectrometry and classical histology on a single tissue section, F. Deutskens, J. Yang, R. M. Caprioli. *Journal of Mass Spectrometry* 2011, 46, 568-571.