

---

# Sublimation/Rehydration of SA or CHCA for imaging proteins/peptides

## Health and Safety:

Wear nitrile gloves, laboratory coat and safety glasses. The process should be done with fume hood ventilation.

## Equipment and Materials:

### Equipment

Petri dish (plastic, Fisher, 08-757-100D, 100mm x 15mm); stainless steel plate (AB, 4347680); heat conductive copper tape (Electron Microscopy Sciences, 77801); Tape (Fisher, 35901R); Oven set at 85°C.

### Chemicals

Acetic acid, methanol; MilliQ water; Carnoy's fluid (6 ethanol, 3 chloroform, 1 acetic acid)

## Procedure

Tissue preparation prior to sublimation.

1. Following tissue sectioning and placement on plate or slide, it may be necessary to wash tissue.
2. Washing conditions for detection of proteins: 70% ethanol (30s), 100% ethanol (30s), Carnoy's fluid (2min), 100% ethanol (30s), H<sub>2</sub>O (30s), 100% ethanol (30s). The section is sublimated with SA to obtain a coating at optimized 0.25mg/cm<sup>2</sup>. (MSRC-R-025)
3. For sections to detect peptides, the section is not washed and is sublimated with CHCA at 0.2 mg/cm<sup>2</sup>.

### Rehydration of sublimated matrices (Figure 1).

1. The slide with sublimated matrix is attached to a stainless-steel plate (used as a heat sink).
2. The plate is attached to the underside of the top part of the Petri dish using a heat conductive copper tape.

3. This top part of Petri dish is placed into a pre-heated oven at 85°C for 2 min.
4. A piece of filter paper is placed in the bottom part of the Petri dish, for SA: 1 mL of MilliQ water and 100 µL of acetic acid or methanol pipetted onto the paper to subsequently create a vapor for the recrystallization process; for peptide analysis using CHCA: 200 µL MilliQ water, 50 µL TFA are pipetted onto the paper.
5. The two parts of Petri dish are reassembled to form a hydration chamber, sealed using a tape (Fisher, 35901R) and left in the oven for 3.5min.
6. The slide is taken out from the petri dish and allow to dry.




Figure 1 Construction of rehydration/recrystallization chamber with a Petri dish, magnet, stainless steel plate and tape.

**Expected outcome:**

The sign of a successful rehydration is determined by a quick color change at the moment the Petri dish is opened. This phenomenon is caused by the evaporation of solvent condensed on the surface of matrix.

**Reference:**

1. Matrix Sublimation/Recrystallization for Imaging Proteins by Mass Spectrometry at High Spatial Resolution, Junhai Yang and Richard M. Caprioli *Analytical Chemistry* **2011** 83 (14), 5728-5734
2. High spatial resolution imaging mass spectrometry and classical histology on a single tissue section, F. Deutskens, J. Yang, R. M. Caprioli. *Journal of Mass Spectrometry* 2011, 46, 568-571.