Algorithm for the Infusion Rate of Glucose During an Insulin Clamp

Jason Blohm, Nicholas Diehl, Joe Jeffrey, Sheng-Yau Lim *GlucoReg*

Background: Hyperinsulinemic clamp

- Plasma insulin rate is raised and held at 120 $\mu U/mL$ via constant insulin infusion
- Glucose is infused and blood glucose levels are measured; GIR adjusted
- Measures insulin sensitivity in patient
- Patients: Pre-diabetics, diabetics, those with endocrine & metabolic disorders



Refresher: Problem Statement

- In hyperinsulinemic clamp studies, Dr. Luther adjusts GIR on the fly based on his clinical judgment
- This can lead to inaccurate adjustments which can affect subject safety and data validity
- Some people claim that an algorithmic approach works, but no one Dr. Luther has talked to has been successful
- We will develop an algorithm that allows researchers to perform these studies in a more controlled manner



Refresher: Needs Assessment \rightarrow Provider

Interface

- 1. Should be simple to understand
- 2. Should include inputs for all possible variables the physician may want to change: target glucose level, insulin clamp level, demographic data, time of experiment
- 3. Given patient demographics and history, should simulate the glucose level over time, prior to clinical testing



Refresher: Needs Assessment \rightarrow Provider

Algorithm

- 1. Should calculate the amount of glucose uptake based on the constant insulin infusion rate the physician specifies
- 2. Depends on actual glucose infusion rate (not suggested)
- 3. Should output a recommended glucose infusion rate that accounts for the time delay in measuring glucose level from blood sample (t-1)



Refresher: Needs Assessment \rightarrow Provider

Timing

- Runtime -- should provide physician with proper glucose infusion rate (GIR) within 10 seconds of inputting the current glucose level
- 2. Should include an easy to navigate UI for immediate data entry



Refresher: Needs Assessment → Patient

Safety

- 1. Ensure that glucose levels do not exceed or drop below safe levels, as determined by the physician
- 2. Measurements need to be taken every 5 minutes to ensure glucose levels are where they should be. If not, the program should alert the physician (future iteration)
- 3. Must run smoothly so that no bugs interrupt the program



Refresher: Needs Assessment \rightarrow System

Applicability and cost

- 1. Should be applicable to different physicians and different hospitals performing the same studies
- 2. Should be open source
- 3. Results from these studies should lower healthcare costs in the future



Progress: Multiple Regression Analysis

Linear regression model:

GIR ~ 1 + Height + Weight + BSA + Age + Gender + Race

Estimated Coefficients:

	Estimate	SE	tStat	pValue
(Intercept)	3.2236	2.011	1.603	0.11325
Height	15.35	5.7597	2.6651	0.0094679
Weight	0.19485	0.083456	2.3348	0.022309
BSA	-20.895	8.0997	-2.5797	0.0119
Age	-0.014226	0.0054871	-2.5927	0.011498
Gender	0.14254	0.16286	0.87521	0.38433
Race	-0.32696	0.22619	-1.4455	0.15261



Progress: First Iteration

- Relying on current algorithm for first iteration of algorithm
 - GIR = ((Gd-Gi)(10)((0.19)(body weight))/((Ginf)(15)) + (((SMi-2)(Gd)(FMi-1))/Gi)
 - Gd = desired plasma glucose concentration (mg/dL), 95
 - Gi = actual plasma glucose concentration (mg/dL)
 - Body weight in kg
 - Ginf = glucose concentration in infusate (mg/mL)
 - SMi = metabolic component: SMi = (SMi-2)(FMi)(FMi-1)
 - FMi = Gd/Gi
 - Will use our linear regression model to set the initial GIR



Progress: GUI

- Will take inputs and generate a predicted GIR
- Once closed, file will save to MATLAB workspace and an excel

file



Code: GUI

% Plot the GIR and GL

```
plot(app.glucoseinfus, time, gir,'-o');
plot(app.glucoseconc, time, gl,'-o');
```

```
app.glucoseinfus.YLim = [0 10];
app.glucoseconc.YLim = [75 105];
```

% Update the time

app.CurrentTime.Value = curtime + 5;

% Update the Infusion Rate

app.InfusionRate.Value = prevgir;

% Button pushed function: Button function ButtonButtonPushed(app, event)

% Update number of button pushes

app.buttonpush.Value = app.buttonpush.Value +1; butpush = app.buttonpush.Value;

% initialize variables

```
studyid = app.StudyID.Value;
height = app.Height.Value;
bsa = app.BSA.Value;
weight = app.Weight.Value;
age = app.Age.Value;
prevgir = app.PrevGIR.Value;
gl one = app.GlucoseLevel1.Value;
gl two = app.GlucoseLevel2.Value;
avg_gl = mean([app.GlucoseLevel1.Value, app.GlucoseLevel2.Value]);
app.AverageGL.Value = avg_gl;
curtime = app.CurrentTime.Value;
LD = app.LowDose.Value;
```

% Calculate the future GIR

```
if prevgir==0 && curgl==0
    app.PredictGIR.Value = 1+ 15.35*height + .19485*weight - 20.895*bsa - age*.014226;
else
    % This is where we will put in the algorithm
    app.PredictGIR.Value = 0;
end
% Update table
app.UITable.Data{butpush,1} = curtime;
app.UITable.Data{butpush,2} = avg_gl;
app.UITable.Data{butpush,3} = prevgir;
if LD == 1
    app.UITable.Data{butpush,4} = 'Low Dose';
else
    app.UITable.Data{butpush,4} = 'High Dose';
end
```

Progress: Exporting to Excel

- Currently, data is manually recorded and transferred to Excel Sheet Template
- Implement code into GUI that exports data when closing out of GUI
- Add button to export to Excel without closing out of GUI



Gantt Chart

Glucose Regulation Algorit...

Preliminary Brainstorming Familiarization and meeting Dr. Luther Brainstorming **Preliminary Data Analysis** Data storage and sorting Database and demographic analysis Multiple regression analysis Algorithm - first round Development of algorithm - first draft Clinical testing (1) Algorithm iterations (1) Algorithm - second round Implementation of machine learning Clinical testing (2) Algorithm iterations (2) **Pump interface** Familiarization with pump/Brainstorm... Prototype First round of testing Prototype modification Finalization

Algorithm + Pump Finalization Creation of user interface Finalization of algorithm, UI, and pum... Design day preparation



Our pump

CareFusion BD Alaris Pump Module with Guardrails MX software suite

Offset Offset Million
August and august
Application Manage (MAPP) Construct Advance Manage (MAPP) Manage
Annual (mark) Annual (
Number (NAME) Auge
International analysis Analysis Like Like <thlike< th=""> Like Like <thl< td=""></thl<></thlike<>
International Sources Internation Sources Internation Sources Internaternation Sources <thinternaternation so<="" td=""></thinternaternation>
Constant Early L00 L00 <thl00< th=""> L00 L00 <thl0< td=""></thl0<></thl00<>
International Configure Date Dir Dir <thdir< th=""> Dir <thdir< th=""></thdir<></thdir<>
Contract Carryon Loss Loss H H H Contract Georgia Loss <
Open Date by Law Table Date Date <thdate< th=""> <thdate< th=""> <thdate< th=""></thdate<></thdate<></thdate<>
Allowing Lass Lass <thlass< thr=""> Lass Lass</thlass<>
Automatic Contract Co
March Marting Average Marting
And Added
Constants Constants 1 Ann
And Lat Lat Lat Lat Lat
Destant Bartists Libs 540 Million
17 84 M 14 14
1



Pump Functionality

- Central side of controller will allow us to add units (sets) for infusion
- Ability to infuse four different drugs at different rates
- Can add four different infusion units for individual infusions
- Can continuously or intermittently deliver infusions
- Guardrails software helps to reduce risk of medication error
- Interface between pump module and PC unit



Pump Concerns (Dr. Luther)

- How to make sure that Dr. Luther is controlling just the one pump and not a pump down the hall
- Getting check-off or approval from VUMC device personnel (will need to get them involved early)
- Once the connection between CPU and pump is made, there needs to be a requirement for manual override should something go wrong
- FDA clearance or approval for testing phase



Pump interfacing: External vs. internal

- Directly with Alaris pump (external)
 - Microcontroller controlled by an external circuit that interfaces with pump tubing coming from Alaris
 - Seems to be repetitive (since Alaris system controls pump as is)
- Internally programming Alaris pump itself
 - If able to take advantage of software, could potentially implement our algorithm within the pump itself



Pump interfacing: Start new

- Microcontroller that controls glucose source directly
 - Syringe pump
 - Potentially program using Arduino





Next Steps

- Familiarization with machine learning
 - Meeting with Dr. Kunda a week from next Wednesday
- Finalize first iteration of algorithm and use in clinical trial (not directly on patient, but side by side)
- Obtain new sources of data
- Begin working on pump interface design/prototype (meeting with Dr. Luther after this)



Questions?

